

Spectra & Filter Efficiency Calculation

Inhaltsverzeichnis

1 The main program - Spectra.....	2
1.1 Spectral Data Files.....	3
2 Main Functions.....	3
2.1 File.....	3
2.2 Dyes.....	3
2.3 Filter.....	4
2.4 Light Source.....	4
2.5 Miscellaneous Spectra.....	4
2.6 Fura Excitation.....	5
2.7 Laser Lines.....	5
2.8 Raman Lines.....	5
2.9 Filter Efficiency.....	5
2.9.1 Scaling Options.....	6

1 The main program - Spectra


In order to use the software one must install Python and some additional modules. The easiest way is probably to install and complete Python distribution called PythonXY. It can be found here:

<http://www.pythonxy.com/>

Actually required modules are:

- Python 2.6
- numpy
- scipy
- wxPython
- matplotlib

Once everything is installed the program can be started just by double-clicking the python script called *spectra.py* and you should see this:



1.1 Spectral Data Files

All used spectral data files must have the following specs:

- 1st column → wavelength [nm]
- 2nd column → transmission / intensity
- decimal separator = . (point, not comma)

All those files are actually text files. The file extensions were adapted:

- *.abs → dye absorption
- *.flu → dye emission
- *.ex → excitation filter
- *.di → dichroic mirror
- *.em → emission filter
- *.ls → light source
- *.txt → any spectral file

2 Main Functions

The main purpose of spectra is to easily display spectral data of dyes and filters to enable the user to select the appropriate “dye/filter/light source” combination.

2.1 File


- *Save Graph As* or *Save Efficiency Graph As* → opens save dialog
- *Set Working Directory* → allows to set a directory (where the dyes & filter are)
- *Quit* → ends the program

2.2 Dyes

- *Load Dye Absorption / Emission* → displays *.abs or *.flu spectral data file
- *Delete Dye Absorption / Emission* → deletes absorption or emission data
- works for 2 sets of dyes

2.3 Filter

- *Load Excitation / Dichroic Mirror / Emission Filter* → displays *.ex or *.di or *.em spectral data file
- *Delete Excitation / Dichroic Mirror / Emission Filter* → deletes the data
- works for 2 sets of filters
- allows simulation of 1 bandpass filter



2.4 Light Source

- *Load Light Source* → displays *.ls spectral data file
- *Delete Light Source* → deletes light source from graph page

2.5 Miscellaneous Spectra


- *Load Miscellaneous Spectra* → displays *.txt spectral data file
- *Delete Light Source* → deletes miscellaneous spectral data from graph page

2.6 Fura Excitation

- displays Fura 340/380 excitation bar to check for filter compatibility

2.7 Laser Lines

- displays various laser lines with the corresponding color




2.8 Raman Lines


- displays possible Raman laser lines for glas and water

2.9 Filter Efficiency

- allows an „rough“ estimation of the detection efficiency
- a separate dialog windows allows choosing the required spectral data



- to see the calculated curves the *Efficiency Graph* must be in focus
 - fluorescence emission spectra is filtered by dichroic and emission filter
 - the areas under the curves are compared



2.9.1 Scaling Options

- Allows adjusting the left and right border of the *Graph* and *Efficiency Graph*